



Neo-Vista

For easy identification of veins in the Newborn



The Neo Vista is a portable vein illumination device for neonates and babies.

When activated it displays a map of the vasculature on the surface of the skin, allowing clinicians to verify patency and avoid valves or bifurcations.

It is a more effective venepuncture procedure, with less patient discomfort, resulting in higher patient satisfaction and ensuring shorter procedure time.

- Easy to use, economical and hands free
- Increases the success rate with first attempt of venepuncture
- Splash proof, so that the user can clean it with disinfectant
- Minimum maintenance required

SIZES AVAILABLE

When ordering please quote the following (minimum order qty 1 box):

ORDER#	DIMENSIONS	QTY/box
TRANS1	108mm x 32mm (max) x 38mm (max)	10

Readmed

Room 2706
Concens Science & Tech Plaza
No. 2911 N. Zhongshan Road
Putou District
Shanghai 200063. P.R.C
Tel: (+86) 21 5013 5443
Email: contact@readmedgroup.com

Usage Instructions

Important: Please read all the instructions, warnings and information below carefully, ensuring that you fully understand them prior to using the Neo-Vista transilluminator.

1. Insert the batteries into the compartment of the Neo-Vista.
2. To activate, press the On/Off button.
3. Hold the Neo-Vista or put on a flat surface and place the palm of the patient's hand over the LED of the device.
4. Dimming the room lights will aid in detection using the Neo-Vista.
5. Insert the IV-Catheter then switch off the Neo-Vista.

N.B. The patient's fingers can be draped over the Neo-Vista in order to stretch the skin, so assisting with the location of veins.

Information:

1. Replace old batteries with new AA type batteries once exhausted, it is recommended that both batteries are replaced at the same time.
2. To clean the Neo-Vista, please refer to the relevant hospital's regulations and protocols that apply.
3. Once the Neo-Vista is no longer operational i.e. LED display faulty and/or On/Off switch broken, ensure disposal is achieved according to the hospital's guidelines.

General Specifications:

Battery type:

AA Alkaline

Compliance:

Standard: IEC 60601-1 General requirements for Medical Equipment

LED type:

Ø5 CAP LED

Output Light wavelength: 620 nm - 635 nm

Technical Specifications:

Lighting Source	Ø5 CAP LED
Working Mode	On-Off
Lumen	7LM
Runtime	12h
Peak Beam Intensity	26 cd
Beam Distance	5m
Impact Resistance	1m
Battery Type	2 x AA Alkaline
Switch Model	Push Button
Output Power	0.12W
Batteries removable	Yes
Dimensions (mm)	10.8 x 3.8 x 3.2cm
Net Weight	26.5g
Material construction	ABS Plastic
Finishing	abrasive
Lens Dimension (mm)	24 x 15.8 x 6.3mm
Lens Material	PC
Switch Position	Body
Logo Position	Body and battery cover

