

Spo2 sensor

TRIMUS KFT.

E-mail : trimus@t-online.hu
 website : www.trimus.hu
 TEL : [06/20-9838-646](tel:06209838646) , Fax: [06/1-410-1148](tel:0614101148)

New Disposable Spo2 sensor (Reusable Spo2 cable + Disposable tape)

BCI

P/N	Description(Connector, Length)
RD001-3	BCI New Reusable cable (with diodes) + Disposable Adhesive tapes sets ,7pin male redel , L=3M



Suitable for BCI 3100, 3101,6100, 9100 Advisor ,
 Autocorr, MiniTorr, Minitorr Plus (**Direct connect**)

P/N	Description(Connector, Length)
RD001-1D	BCI DB9 Reusable cable (with diodes) + Disposable Adhesive tapes , L=1M

Compatibilities: Compatible with BCI 3300, 3301, 3302,3303, 3304, 3401



Biosys (Korea)

P/N	Description(Connector, Length)
RD002	Biosys New Reusable cable (with diodes) + Disposable Adhesive tapes sets ,MD6M , L=3M

compatible with **Biosys** BPM-300



CSI

P/N	Description(Connector, Length)
RD004-3	CSI Reusable cable (with diodes) + Disposable Adhesive tapes sets , 5pin Lemo male L=3M

:Criticare Systems International :CSI 501, 503, 504, 507, 508, 581,1100, Poet (602-3A),8100

PONI,Scholar (**Direct connect**)




CSI

P/N	Description(Connector, Length)
RD005-1	CSI Reusable cable (with diodes) + Disposable Adhesive tapes sets , DB9 r

Compatibilities :Criticare Systems International :CSI 503spot,504plus,504DX, 504DXN.


(**Direct connect or use Adapter EX005**)

Spo2 sensor



Critikon

P/N	Description(Connector, Length)
RD006	CSI Reusable cable (with diodes) + Disposable Adhesive tapes sets , 7pin Lemo male L=3M



Compatibilities : Suitable for Critikon Dinamap 8700, 8710, 8720 (direct connector)

Datex -Ohmeda

P/N	Description(Connector, Length)
RD007-1	Datex Reusable cable (with diodes) + Disposable Adhesive tapes sets , 7pin Lemo male L=3M



Datex-Ohmeda

P/N	Description(Connector, Length)
RD007-3	Datex Reusable cable (with diodes) + Disposable Adhesive tapes sets , 10pin female, L=3M

Compatibilities: Suitable for **Datex-Ohmeda** As/3,CS/3,CH/S,CH/2S,CH-RS,CG/2S,Cardiicap II , Capnomac Ultima, light monitor , m-Estpr. Oscar II, OSP-200
Satlite Trans, Satlite Trans Plus , ULTSVI, ULTSVO,ULTSIO,ULT200.(Direct Connect)



Datascope

P/N	Description(Connector, Length)
RD008-3	Datascope Reusable cable (with diodes) + Disposable Adhesive tapes,8pin male,L=3M



Compatibilities : Datascope Accusat ,
Accutorr Plus , passport, Accutorr 3/4 SAT,
(Direct connect)

Datex Ohmeda

P/N	Description(Connector, Length)
RD009	Ohmeda Reusable cable (with diodes) + Disposable Adhesive tapes,hyperstronic 7pin ,L=3M

Compatibilities: Be Compatible with **Ohmeda** 3700,,3710,3770,3775,3800,3900,Biox3740, 3760, 4700, Oxicap
GE Medical Systems :7000 series ,Corometricsa 511,Tram A.R.,TRAM Modular X00sl series
Mennen :Horizon XL (Direct connect)



Spo2 sensor

Datex-Ohmeda

P/N	Description(Connector, Length)
RD051	Datex-Ohmeda Reusable cable (with diodes) + Disposable Adhesive tapes,8pin male ,L=1M .

Compatibilities : Use with Handheld pulse oximeters (Tuffsat, 3775)
 Stand Alone pulse oximeter (37** series, 3800, 3900.3900P),respiratory gas monitors (4700, 5259RGM)
 M-OSAT module , N-XOSAT option for Cardiocap.5 ; OEM patnering monitors ;
 All S/5 monitors, all AS/3,CS/3 Monitors , Cardiocap ,Capnomac Utima and Satlite
 Require Adapter EX043 or EX044 or EX042



Datex-Ohmeda

P/N	Description(Connector, Length)
RD063	Datex-Ohmeda Trusat Reusable cable (with diodes) + Disposable Adhesive tapes,DB9 male,L=1M .



GE-Marquette

P/N	Description(Connector, Length)
RD039	Datex-Ohmeda Reusable cable (with diodes) + Disposable Adhesive tapes,11pin male,L=3M .


Suitable for **GE Medical Systems** : Dash,Ge Medical System : Solar, eagle
 GE Medical System : TRAM 250,(Direct Connect)



GE-Ohmeda

P/N	Description(Connector, Length)
RD056	Datex-Ohmeda Reusable cable (with diodes) + Disposable Adhesive tapes,11pin male,L=3M .

Compatible with Datex-Ohmeda S/5 FM Port ,GE Anesthesia monitor



INVIVO

P/N	Description(Connector, Length)
RD012	Invivo Reusable cable (with diodes) + Disposable Adhesive tapes, Lemo 7pin male,L=3M .


Compatibilities:Use with **Invivo** : 4500 Scout,
 Scout+/-1/+2/+3/+4 , Scout es , Millennia, 4500 Plt
 Omega 1445, Omni-Trak TVS, Sat Scout
Direct Connect .

Spo2 sensor



MEK


P/N	Description(Connector, Length)
RD013-1	MEK Reusable cable (with diodes) + Disposable Adhesive tapes, DB9 male,L=3M .



Compatibilities :Use with Mek patient monitor
MP-100,MP-110, MP-400,MP-1000

Novamatrix

P/N	Description(Connector, Length)
RD023	Novamatrix Reusable cable(with diodes)+Disposable Adhesive tapes, Hyperstronics 7pin ,L=3M .



Compatibilities :Suitable for **Novamatrix**

505,510,511,515B/C,520, 7000, 7100, Oxypleth

Spacelab Medical : 90351-D/G, 90465,9.466,90467,
90489,90651-07/0890419-01,90651A-07 (**Direct Connect**)

P/N	Description(Connector, Length)
RD060	Novamatrix Reusable cable (with diodes) + Disposable Adhesive tapes, DB9 male ,L=3M .



Nellcor

P/N	Description(Connector, Length)
RD024-1	Nellcor Reusable cable (with diodes) + Disposable Adhesive tapes,DB9 male ,L=1M .
RD024-3	Nellcor Reusable cable (with diodes) + Disposable Adhesive tapes,DB9 male ,L=3M .



Compatibilities : Suitable for (**Artema** : Diascope 2, Gamme 8000, **Biosys** BPM200, **Drager** vitalert 1000, Cicero Narkomed, Vitara,PM8014, D8060, narkomed, **Ge-Marquette** dashEagle 3000, 4000, solar, TRAM X50-Modules, **Hellige** Eagle 1000,2000,3000 **HP** 78352A/C, 78254A/C,Merlin:1020A,1024A,1025A/B 1205A,1722A/B, 1732A/B,2475B, **Invivo Research** 4500 scout, 4500 plus, **Ivy Biomedical** 405A, 405D, **Kontron** Micromon, Supermon ,Kolormon , 7840,7845, **Pach tech** vitalmax, 800series Schiller CM8, **Nellcor** N100, N200, N180, N6000, N190, N20, N3000. Require Adapter cable , **Datascope** : Expert , Passport **Require Adapter EX020 or EX022 or EX021**

Nellcor OXIMAX

P/N	Description(Connector, Length)
RD024-1-OXI	Nellcor Reusable cable (with diodes) + Disposable Adhesive tapes, DB9 with

Created with


Spo2 sensor

Compatibilities Compatible with **GE-MEDICAL SYSTEM**: DINASMAP, Pro1000;
NELLCOR : N550, N595 . Require Spo2 Adapter Cable EX023 .



Nonin

P/N	Description(Connector, Length)
RD032-1	Nonin Reusable cable (with diodes) + Disposable Adhesive tapes, DB9 ,L=1M .



Compatible with Nonin 8500, 8600, 8700, 8800.

P/N	Description(Connector, Length)
RD036	Nonin Reusable cable (with diodes) + Disposable Adhesive tapes,8pin female,L=3M .

Compatible with Nonin 8600 Pulse oximeter



P/N	Description(Connector, Length)
RD058	Nonin Reusable cable (with diodes) + Disposable Adhesive tapes,Lemo 6pin male ,L=3M .



Compatible with **Nonin 8600 Pulse oximeter**

Nihon Kohden

P/N	Description(Connector, Length)
RD022-3	Nihon Kohden Reusable cable (with diodes) + Disposable Adhesive tapes,10pin round ,L=3M .



Compatibilities : Use with ***Nihon Kohden*** :
 AL-800PA,Module,BSM-1101, 2101/2,3101,7100/7200
 OLV-1100A,lifescop , BSM-8200/1 Lifescop 12 ,
 BSM-8301/2 Lifescop 9,SM-8800 Lifescop 14

P/N	Description(Connector, Length)
RD022-1	Nihon Kohden Reusable cable (with diodes) + Disposable Adhesive tapes, DB9 male ,L=1M .

Compatibilities : Use with ***Nihon Kohden*** :BSM-2300 Lifescop ; BSM-4100 Lifescop ;
 BSM-5100 Lifescop ;BSM-9510 Lifescop , BSM-9800 Lifescop
 Require Spo2 Adapter cable EX041



Spo2 sensor

Nihon Kohden

P/N	Description(Connector, Length)
RD052	Nihon Kohden Reusable cable (with diodes) + Disposable Adhesive tapes,14pin male ,L=3M .

Use with **Nihon Kohden** :BSM-2300 Lifescope i/L; BSM-4100 Lifescope P;
BSM-5100 Lifescope A ;BSM-9510 Lifescope M , BSM-9800 Lifescope s
OPV-1500 Lifescope N, WEP-4200, WEP-8400 , DDG-2001



Philips /HP

P/N	Description(Connector, Length)
RD025	HP Reusable cable (with diodes) + Disposable Adhesive tapes, Round 12pin ,L=3M .

Compatibilities: Suitable for Philips 78352A/C,78354A/C,78834C

Miracle system, Viridia A1/A3 (**Direct Connect**)



P/N	Description(Connector, Length)
RD026	HP Reusable cable (with diodes) + Disposable Adhesive tapes, 8pin male,L=3M .

Compatibilities :Use with **HP** ,Viridia ,M3, M4, M2601A
M3000A, M3500B, M4735A (**Direct connect or require adapter EX025, EX026**)



Palco

P/N	Description(Connector, Length)
RD027-3	PALCO Reusable cable (with diodes) + Disposable Adhesive tapes,RJ12-6P6C,L=3M .

Compatibilities Suitable for Palco Monitor (Direct Connect)



S&W artema

P/N	Description(Connector, Length)
RD029	S&W Artema Reusable cable (with diodes) + Disposable Adhesive tapes, 10pin female ,L=3M .

Compatibilities:Suitable for S&W artema (diascope G2)Patient monitor (Direct Connect)



P/N	Description(Connector, Length)
AF073	S&W Artema Reusable cable (with diodes) + Disposable Adhesive tapes, 12c

Spo2 sensor



Pace Tech

P/N	Description(Connector, Length)
RD028	Pace Tech Reusable cable (with diodes) + Disposable Adhesive tapes, Switchcraft 5pin ,L=3M .

Suitable for Pace Tech all with / Black Face (Except Vitalmax 88 w/ Nellcor)



Mindray

P/N	Description(Connector, Length)
RD016	Mindray Reusable cable (with diodes) + Disposable Adhesive tapes, 6pin Redel male (dual), L=3M .




Compatible with **Mindray** 9000 (**Direct Connect**)

P/N	Description(Connector, Length)
RD017	Mindray Reusable cable (with diodes) + Disposable Adhesive tapes, 5pin Redel (single), L=3M .

Compatible with **Mindray** PM9000, **Goldway** UT4000

Pace Tech , **Biolight** , **Anke** , **ChoiceMed** .



P/N	Description(Connector, Length)
RD018	Mindray Reusable cable (with diodes) + Disposable Adhesive tapes, 6pin REDEL (dual) , L=3M .



Suitable for **Mindray** PM9000 monitor (**Direct Connect**)

Also compatible with Creative PC9000

P/N	Description(Connector, Length)
RD019	Mindray Reusable cable (with diodes) + Disposable Adhesive tapes, 6pin REDEL (dual) , L=3M .



Suitable for **Mindray** PM8000 monitor

Also Compatible with Edan (Shenzhen)

P/N	Description(Connector, Length)
RD020	Mindray Reusable cable (with diodes) + Disposable Adhesive tapes, lemo 5pin, L=3M .

Suitable for **Mindray** PM9000/509B/500A monitor (Direct Connect)



P/N	Description(Connector, Length)
RD021	Mindray Reusable cable (with diodes) + Disposable Adhesive tapes, 6pin REI

Spo2 sensor



Suitable for Mindray PM600 (**Direct Connect**)

P/N	Description(Connector, Length)
RD042	Mindray Reusable cable (with diodes) + Disposable Adhesive tapes, 12pin male ,L=3M .



P/N	Description(Connector, Length)
AF018 -M	Mindray cable(with diodes) + Disposable Adhesive tapes,,6PIN REDEL Male(dual,60 degree), L=3M.

**PM-7000, PM-8000, PM-9000, MEC-1000, PM-9000Express,
PM-600, new PM-6000, VS-800, (After October ,2005)**

With Masimo Technology



P/N	Description(Connector, Length)
RD077	Mindray Reusable cable(with diodes) + Disposable Adhesive tapes ,Round 7pin, L=3M.

Be Compatible with Mindray Beneview T5/T6/T8

With Nellcor Technology



Spacelab :

P/N	Description(Connector, Length)
RD030	Spacelab Reusable cable (with diodes) + Disposable Adhesive tapes,Hyperstronics 12pin ,L=3M .



Compatibilities: Suitable for Spacelab :90419-01,
90651A-07,90351-D/Gm; 90465,90466,90467,90489,
90501 (**Direct connect**)

New Spacelab


P/N	Description(Connector, Length)
RD031	Spacelab Reusable cable (with diodes) + Disposable Adhesive tapes, Hyperstronics 12pin ,L=3M .



Suitable for Spacelab :New modules 90651-08, 90496, Ultraview (Direct Connect)

Siemens

P/N	Description(Connector, Length)
RD040	Siemens Reusable cable (with diodes) + Disposable Adhesive tapes,10pin male ,L=3M .



Compatibilities : Suitable for Siemens Medical:
Sirecust 700, Siemens Medical: Siresust 630 , Sirem Spo2 module

Spo2 sensor

P/N	Description(Connector, Length)
RD059	Siemens Reusable cable (with diodes) + Disposable Adhesive tapes, Siemens 7pin ,L=3M .



Compatibilities :Suitable for **Siemens Medical** :
SC 6002XL, SC7000, SC8000, SC9000XL

Dolphin

P/N	Description(Connector, Length)
RD041	Dolphin Reusable cable (with diodes) + Disposable Adhesive tapes, DB9 male,L=1M .

Compatible with Dolphin Pulse oximeter 2150, 2100

LOHMEIER

P/N	Description(Connector, Length)
RD048	LOHMEIER Reusable cable (with diodes) + Disposable Adhesive tapes,10pin ,L=3M .



Compatible with LOHMEIER M010 .

Simed /Baxter

P/N	Description(Connector, Length)
RD049	Simed Reusable cable (with diodes) + Disposable Adhesive tapes, Hyperstronics 9pin ,L=3M .



Be Compatible with **Simed** S100, S50, S100E

Corometrics

P/N	Description(Connector, Length)
RD050	Corometrics Reusable cable (with diodes) + Disposable Adhesive tapes,Hyperstronics 9pin,L=3M .




Compatible with **Corometrics** 118 , 129

Kontron

P/N	Description(Connector, Length)
RD053	<u>Kontron</u> Reusable cable (with diodes) + Disposable Adhesive tapes, Round 6pin ,L=3M .

Comptibilities

Use with Kontron Micromon, Minimon, Supermon, Kolormon, 7840, 7845,
7137




Nellcor/Colin (D)

P/N	Description(Connector, Length)
RD054	Nellcor / Colin Reusable cable (with diodes) + Disposable Adhesive tapes,14

Created with

Spo2 sensor



Compatibilities : Suitable for **Datascope** : Accutorr Plus,
Nellcor : NPB-290/5, N-390, N-395, N-3000,
 N-6000, NPB-290-5, Colin BP88/306

Goldway

P/N	Description(Connector, Length)
RD055	Goldway Reusable cable (with diodes) + Disposable Adhesive tapes, lemo 5pin ,L=3M .



AF055 is compatible for **Goldway**: 4000F
Direct Connect (直接接机器端)

AF017 (Mindray Compatible) ALSO COULD BE USE FOR GOLDWAY UT4000A

Schiller

P/N	Description(Connector, Length)
RD064	Schiller Reusable cable (with diodes) + Disposable Adhesive tapes, 7pin ,L=3M .

Compatible with Schiller Argus Pro . Argus LCM

M&B Electronic

P/N	Description(Connector, Length)
RD066	M&B Reusable cable (with diodes) + Disposable Adhesive tapes, Hyperstronics 9pin ,L=3M .

Compatible with M&B : CD2000 , MB526T



Comen

P/N	Description(Connector, Length)
RD065	M&B Reusable cable (with diodes) + Disposable Adhesive tapes, Redel 6pin (Dual) ,L=3M .



Use with Comen Star-1000

FMT

P/N	Description(Connector, Length)
RD068	FMT Reusable cable (with diodes) + Disposable Adhesive tapes, Redel 6pin (Dual) ,L=3M .



Bruker

P/N	Description(Connector, Length)
AT003	Bruker Reusable cable (with diodes) + Disposable Adhesive tapes, Round 12pin ,L=3M .

Compatibilities : **Compatible with Bruker (Odam): All**



Spo2 sensor

Mennen

P/N	Description(Connector, Length)
RD062	Mennen Reusable cable (with diodes) + Disposable Adhesive tapes, Round 13pin ,L=3M .



Compatible with Mennen Monitor Envoy
VITALOGIK 5000/ 5500 Model

P/N	Description(Connector, Length)
RD072	LNCS DC-1 Reusable cable (with diodes) + Disposable Adhesive tapes, DB9 ,L=3M .



OEM part No. : LNCS DC-1
Connect with extension cable LNC-04, LNC-10, LNC-14
0515-30-11221 use with Mindray Old PM6000
0010-30-42625 use with Mindray PM9000
LNC-10-GE (2027263-002) use with GE medical
PM-6000series, PM-7000, PM-8000series, PM-9000series
VS-800, BeneView

Neusoft

P/N	Description(Connector, Length)
RD075	Neusoft Reusable cable (with diodes) + Disposable Adhesive tapes,Redel 6pin ,L=3M .



P/N	Description(Connector, Length)
RD075	Newtech Reusable cable (with diodes) + Disposable Adhesive tapes,8pin ,L=3M .



Use with Newtech NT2A

P/N	Description(Connector, Length)
RD081	General Meditech cable (with diodes) + Disposable Adhesive tapes, DB9 male, L=1M



With digital Technology
Be compatible with General Meditech G1A



KEMP

P/N	Description(Connector, Length)
RD082	KEMP Reusable cable (with diodes) + Disposable Adhesive tapes, Redel plastic 6pin , one bit



Use with KEMP: KP-900C/900B

Spo2 sensor

GuoTeng/ PuAo

P/N	Description(Connector, Length)
RD083	GuoTeng Reusable cable (with diodes) + Disposable Adhesive tapes, Redel plastic 6pin , one bit



Greatmade P/N	Description(Connector, Length)
DPT-A	Disposable tape separately .

