## 740 SELECT



### 740SELECT ECG

The 740SELECT monitoring platform with 3-Lead ECG is highly configurable, offering best-in-class technology, software and hardware upgrade options.

# **Ordering Information:** 740SELECT ECG monitor configurations with software and hardware options

Part Number	Parameters Installed	Model Number*	<b>Software Options</b> (factory or field enabled)	
			CO2	Masimo rainbow
01-02- <b>1018</b>	NIBP, Masimo SpO2, ECG	XM0E0		SpCO®, RRa™, PVI®
01-02- <b>1019</b>	NIBP, Nellcor SpO2, ECG	XN0E0		
01-02- <b>1020</b>	NIBP, Masimo SpO2, FILAC, ECG	XMFE0		SpCO®, RRa™, PVI®
01-02- <b>1021</b>	NIBP, Nellcor SpO2, FILAC, ECG	XNFE0		
01-02- <b>1022</b>	NIBP, Masimo SpO2, Exergen, ECG	XMAE0		SpCO®, RRa™, PVI®
01-02- <b>1023</b>	NIBP, Nellcor SpO2, Exergen, ECG	XNAE0	Masimo	
01-02- <b>1024</b>	Masimo SpO2, ECG	0M0E0	or Covidien	SpCO®, RRa™, PVI®
01-02- <b>1025</b>	Nellcor SpO2, ECG	0N0E0		
01-02- <b>1026</b>	Masimo SpO2, FILAC, ECG	0MFE0		SpCO®, RRa™, PVI®
01-02-1027	Nellcor SpO2, FILAC, ECG	0NFE0		
01-02- <b>1028</b>	Masimo SpO2, Exergen, ECG	0MAE0		SpCO®, RRa™, PVI®
01-02-1029	Nellcor SpO2, Exergen, ECG	0NAE0		

The 740SELECT is designed for use in a wide variety of patient care settings:

Emergency department

Procedure rooms
EMS/Pre-hospital
Bedside monitoring
Surgical centers
Alternate care
Dental offices
Sub-acute care
Urgent care center
Transport

- X = CASMED MAXNIBP®
- M= Masimo® rainbow SET®
- N = Nellcor™ OxiMax™
- A = Exergen Temporal Artery TemperatureScanner™
- F = FILAC<sup>™</sup> 3000 Predictive Temperature
- E = 3-Lead ECG



Shown with ECG, SpO<sub>2</sub>, NIBP, Temp and EtCO<sub>2</sub> option



Shown with ECG, Resp, SpO<sub>2</sub>

#### 740SELECT 3-Lead ECG Specifications

- Patient Type: adult / pediatric/ neonate
- HR (15-300 bpm) / resolution 1 bpm
- Impedance respiration rate (RR)
- Alarms & alarm messages
  - Upper and lower HR limit - Upper and lower RR limit
  - Asystole, ventricular fibrillation
  - Lead off, artifact
  - (RR) Out of range (too high)
  - HR and RR auto upper and lower alarm limit selection
- Pacer pulse detection yes, lead I, II or III
- Pacer (spike) display: no (default), yes
- · Waveforms:
  - Leads: I, II (default), or III
  - ECG waveform cascade selection: off (default), on
  - Freeze main screen waveforms yes
  - ECG waveform size: 2.5, 5, 10 (default), 15 or 20 mm/mV  $\,$
  - Sweep speed: 6.26, 12.5 or 25 (default) mm/sec
  - RESP (waveform) enabled: no (default), yes
- Audible pulse tone: off (default), on
- HR source selection: auto (default), ECG, SpO<sub>2</sub>, NIBP
- ECG filter: monitor (0.67-40 Hz) (default) and diagnostic (0.05-40 Hz)
- Defibrillation protection: complies with IEC 60601-2-27
- CAL pulse: inserts 1mV pulse into ECG waveform
- Monitor suitable for use in the presence of electro-surgical (ESU) equipment - yes

#### **ECG Start-Up Kits**

AAMI		IEC		
01-02-0891	Adult/Ped	01-02-0893	Adult/Ped	
01-02-0892	Neonatal	01-02-0894	Neonatal	
Start-Up Kit Contents:				
Adult/Ped Kit: Trunk Cable, Snap Lead Set Neonatal Kit: Trunk Cable, Mini-Clip Lead Set ECG Quick Reference Guide included in all kits				

#### **ECG Accessories**

AAMI		IEC		
	01-02-0920	Trunk Cable, 2m	01-02-0929	Trunk Cable, 2m
	01-02-0923	Snap Lead Set, 1m	01-02-0927	Snap Lead Set, 1m
	01-02-0921	Mini-Clip Lead Set, 1m	01-02-0925	Mini-Clip Lead Set, 1m

Contact Zoe Medical customer support for ordering information at customersupport@zoemedical.com.





#### **740SELECT Monitor Features**

#### **Bright 7" Touchscreen Display**

- Easy to read color coded numerics
- Readily operated with a gloved hand
- Auto dim (night mode)
- · Wide viewing angle, anti-glare coating
- Touch click feature for user feedback

#### **Intuitive Operation**

Alarm Limit Adjustment

x% of the current

measured value

Upper

125%

Lower

80%

- One touch access to most used settings and features
- Fixed touch keys for most frequently used functions

#### **Operating Features**

- Screen lock for added patient safety when needed
- Display up to 3 waveforms
- Alarm log displays recent physiologic alarms and setting changes
- Trends 1 minute average of continuously measured parameters
- Standby mode for enhanced bedside workflow / alarm management
- Li-lon battery enhances battery operation time and charging performance
- All-in-one, continuous monitoring and spot check mode
- Sturdy functional handle for portability, ease of carrying and transport
- Battery status indicator
- Software updates via USB

#### **Parameters**

- 3-Lead ECG, HR, Impedance Respiration (RR)
- SpO<sub>2</sub> / PR Masimo rainbow SET or Nellcor OxiMax
- NIBP / PR CASMED Max NIBP **M**otion **A**rtifact e**X**traction Technology
- Temperature FILAC Predictive or Exergen Temporal Artery Scanner

#### Masimo rainbow Parameters

- Acoustic Respiration, RRa
- Carboxyhemoglobin, SpCO
- Pleth Variability Index, PVI

#### **Software and Hardware Options**

	Modular CO <sub>2</sub> Options		
	Sidestream		Mainstream
Capnography	Covidien Microstream™ Micropod™	Masimo ISA™	Masimo IRMA™

#### Full Line of accessories to meet your workflow needs

- Monitor: roll stand, transport, swivel arm, EMS, carry bag
- Bar code reader: patient and clinician IDN
- Thermal Printer: No batteries or external AC power required
  - Print on alarm, print on save
  - 2 waveforms (selectable)

#### Portable, Lightweight, Rugged Design

- 4.7 lbs, 2 kg (with Li-lon battery)
- Tested and complies with the more rigorous standards called out in ISO 80601-2-61:2011

#### Alarms

- Audio and visual
- Upper and lower alarm limits
- Alarm silence and alarm pause
- · Low battery audio and visual alert
- Secondary audio alert (timed)
- Dedicated message area

#### 740SELECT is available with and without ECG



Exergen Temporal Artery Scanner



FILAC 3000 Predictive Temperature

MAXNIBP is a registered trademark of CAS Medical Systems, Inc. 740SELECT is a trademark of Zoe Medical, Inc. Temporal Scanner is a trademark of Exergen Corporation. Nellcor and OxiMax are trademarks of a Covidien company. Micropod and Microstream are trademarks of a Covidien company registered in the U.S. and foreign countries. FILAC is a trademark of Covidien AG registered in the U.S. and foreign countries. Masimo, rainbow SET, SpCO and PVI are registered trademarks of Masimo Inc. IRMA, ISA and RRa are Masimo, Inc. trademarks. Zoe Medical reserves the right to make changes to this document without notice.